

Smallest but most advanced

CP-770

CHART PROJECTOR



CP-770 Specifications

Charts	33 (type M, type T, type G) 27 (type UK)
Refraction distance	2.9 - 6.1 m
Chart magnification (at 5m)	30x (type M, type T, type G), 25x (type UK)
Chart size (at 5m)	275 mm dia., 330 (W) x 225 (H) : type M, type T, type G 250 mm dia., 250 (W) x 225 (H) : type UK
Light source	3W (LED)
Chart brightness	230 cd/m ²
Masks	Open, Horizontal line, Vertical line, Single letter
Filter	Red / Green
Chart changeover devices	Remote control Control box for NIDEK Smart Refractor for RT-5100 (type M, type UK, type T, type G, type pM) for RT-2100/900 (type M, type UK, type T)
Chart changeover speed	0.03 seconds / 1 frame
Auto-OFF function	5 / 15 minutes / Off selectable
Power supply	AC 100 V to 240 V 50 / 60 Hz
Power consumption	70 VA
Dimensions & Weight	
Main body	177 (W) x 278 (D) x 192 (H) mm / 3.5 kg 6.97 (W) x 10.94 (D) x 7.56 (H) " / 7.72 lbs.
Remote control	66 (W) x 184 (D) x 21 (H) mm / 120 g 2.6 (W) x 7.24 (D) x 0.84 (H) " / 0.27 lbs.
Standard accessories	Remote control, Screen, Polarized glasses, Wood screw (x4), Power cord, Batteries (x2), Chart sheet
Optional accessories	Table stand, Floor stand, Wall bracket, Red / Green glasses, Cable, Screen angle fixing set, Mirror

Key Features

- Maintenance-free light for improved cost efficiency
- Brighter and clearer chart display with white LED
- Optical zoomer allows a variable installation distance
- Sophisticated design



Eye & Health Care
NIDEK CO., LTD.

HEAD OFFICE

34-14 Maehama, Hiroishi
Gamagori, Aichi 443-0038, Japan
Telephone : 81-533-67-6611
Facsimile : 81-533-67-6610
URL : <http://www.nidek.co.jp>

TOKYO OFFICE

(International Div.)
3F Sumitomo Fudosan Hongo Bldg.,
3-22-5 Hongo, Bunkyo-ku, Tokyo,
113-0033 Japan
Telephone : 81-3-5844-2641
Facsimile : 81-3-5844-2642
URL : <http://www.nidek.com>

NIDEK INC.

47651 Westinghouse Drive
Fremont, CA 94539, U.S.A.
Telephone : 1-510-226-5700
 : 1-800-223-9044 (US only)
Facsimile : 1-510-226-5750
URL : <http://www.usa.nidek.com>

NIDEK SOCIÉTÉ ANONYME

Europarc
13, rue Auguste Perret
94042 Créteil, France
Telephone : 33-1-49 80 97 97
Facsimile : 33-1-49 80 32 08
URL : <http://www.nidek.fr>

NIDEK TECHNOLOGIES SRL

Via dell'Artigianato, 6 / A
35020 Albignasego (Padova), Italy
Telephone : 39 049 8629200 / 8626399
Facsimile : 39 049 8626824
URL : <http://www.nidektechnologies.it>

*Specifications and design are subject to change without notice for improvement.



Printed on environment-friendly recycled paper.