

TL-228M / 228P

Trial lens sets



MAGNON

TL-228M / 228P

MAGNON

Trial lens sets

The MAGNON TL-228M and TL-228P are double curve trial lens sets which offer a standard selection of lens powers to meet practitioner's need.

FEATURES

TL-228M

- ★ Double curve lenses of 38mm in diameter
- ★ Metal rims (Gold-plated for cylinder lenses)
- ★ Removable plastic inner tray with space for the trial frame
- ★ Leatherette carrying case with lock

TL-228P

- ★ Double curve lenses of 22mm in diameter
- ★ Color coded plastic rims of 38mm in diameter
- ★ Removable plastic inner tray with space for the trial frame
- ★ Leatherette carrying case with lock

SPECIFICATIONS

Spheres	35 pairs each of plus and minus
Cylinders	17 pairs each of plus and minus
Prisms	10 pieces
Accessories	10 pieces
	<ul style="list-style-type: none"> ● : Occluder1 ○ : Pin hole (1mm)1 ⊖ : Slit (1mm)1 ⊕ : Test mark for PD2 R : Red lens1
	<ul style="list-style-type: none"> G : Green lens1 F : Frosted glass1 N : 0.0D lens1 M : Maddox1

- Spheres	- Cylinders	+ Cylinders	+ Spheres
20	6.00	6.00	20
18	5.00	5.00	18
16	4.50	4.50	16
15	4.00	4.00	15
14	3.50	3.50	14
13	3.00	3.00	13
12	2.50	2.50	12
11	2.25	2.25	11
10	2.00	2.00	10
9.00	1.75	1.75	9.00
8.00	1.50	1.50	8.00
7.50	1.25	1.25	7.50
7.00	1.00	1.00	7.00
6.50	0.75	0.75	6.50
6.00	0.50	0.50	6.00
5.50	0.25	0.25	5.50
5.00	0.12	0.12	5.00
4.50			4.50
4.00			4.00
	Accessories	Prisms	
3.75	● ○ ⊖ ⊕ ⊕ R G F N M	0.5° 1° 2° 3° 4° 5° 6° 7° 8° 10°	3.75
3.50			3.50
3.25			3.25
3.00			3.00
2.75			2.75
2.50			2.50
2.25			2.25
2.00			2.00
1.75			1.75
1.50			1.50
1.25			1.25
1.00			1.00
0.75			0.75
0.50			0.50
0.25			0.25
0.12			0.12

Specifications and design are subject to change for improvement without notice.

H. OGINO & CO., LTD.

Yokohama Port P.O.Box 170, Yokohama 231-8691, Japan
Tel: (045)662-5261 Fax: (045)651-5467