

U-1500

Ophthalmic instrument positioning unit

SPECIFICATIONS

CHAIR

Elevation stroke	240 mm (S type) 270 mm (L type)
Seat height	520~760 mm (S type) 550~820 mm (L type)
Range of reclining back	100°~170°
Range of rotation	360°
Weight	105 kg
Power source	AC 100, 120, 220 or 240 V
Power consumption	400 VA
Finish	Blue

STAND

RT arm	Vertical range : 200 mm Weight capacity : 8 kg Diameter of shaft : ϕ 20.5 mm
CP arm (option)	Weight capacity : 8 kg
Sliding table	Slide stroke : 360 mm Power supply : Selectable primary and secondary voltages Weight capacity : 50 kg Height adjustment : 30 or 60 mm collar (option)
Rear panel	Two auxiliary electric sockets, two fuses (10 A)
Weight	119 kg
Power source	AC 100, 120, 220 or 240 V
Finish	White grey and blue
Accessories	Five spare fuses, extension cable for head lamp, two 3-pin connectors and tools.
Optional accessories	Foot switch, CP arm, table height extension collar (30 or 60 mm), support column for HAAG type SL cord, third arm for another instrument or LCD monitor.

*Specifications and design are subject to change for improvement without notice.



H.OGINO & CO.,LTD.

Yokohama Port P. O.Box 170, Yokohama 231-8691, Japan
Tel: (045)662-5261 Fax: (045)651-5467
E-mail : magnon@h-ogino.co.jp URL : www.h-ogino.co.jp

MAGNON