



DRY EYE TREATMENT

CLASS II - HOME HEALTH REMEDY



December 2021
ver. 1 - 2021



OPHTHALMOLOGIST

MDs can use any kind of Dry Eye treatment medical devices such as Lipiflow, IPL or QMR. All of them are very expensive for both patient and doctor.

Activa is the unique, safe, and Effective treatment available for a fraction of the cost of other Dry Eye medical devices.

After each use the mask can be easily sanitized and does not require any gel during treatment as many other devices do.

CLASS II - HOME HEALTH REMEDY

Activa is the first dry eye medical device designed for both in-clinic and home use!


It is clinically proven to improve tear film quality and patient comfort. The Doctor can **customize number of sessions available** for the rental period of your customer to a limited amount according to the prescribed therapy.

A natural way to significantly improve tear stability is to increase Meibomian gland output - this means your eyes stay lubricated and symptom-free longer, as nature intended.

A clinically proven mechanical method to de-cap and stimulate Meibomian glands is applying heat at a very specific temperature and to massage the melted oil out of the glands. This procedure unclogs your glands and stimulate the natural production of oils to keep your eyes lubricated.

Activa is the innovative eye-mask providing a complete heating and automatic squeezing technology in a fast and painless single device.

Patients typically notice a distinct improvement in comfort and a significant reduction in their dry eye symptoms in the first month of use. It is effective on meibomian glands disfunction associated with tear film evaporation.



ARE YOU USING AN IPL IN YOUR CLINIC?



As you may know, **IPL patients are not feeling the effect of the treatment from the first application and meibomian glands expression and squeezing is an important practice combined to the IPL therapy.**

Including Activa in your Evaporative Dry Eye therapy plan combined to the IPL effects will solve the above issues by mechanically warming and massaging the lids of the patient. This will express the meibum out of the glands increasing the efficiency of IPL effects and providing to the patient an immediate feeling of relief and improvement of the pathology.

DON'T YOU HAVE AN IPL IN YOUR PRACTICE?

Activa is a registered medical device, perfect for an efficient treatment of Evaporative Dry Eye disease and offered for one of the lowest prices compared to other available medical devices for dry eye therapy. This makes it the perfect device to start treating your patients with.

- Touch screen controlled, user-friendly and easy to manage.
- Through a fully automated procedure the device is able to melt the meibum inside the glands at 42°C and simultaneously squeeze them, this way releasing the melted oil and allowing a strong improvement of the lipid layer quality.
- The complete procedure takes only 15 minutes and can be performed in any MD/OD practice.



HOW DOES IT WORK

This innovative Eye-mask provides a complete heating and automatic squeezing technology in a single device fast and painless. Effective on Meibomian glands dysfunction associated with tear film evaporation.

- **Unclog glands**
- **Melting the meibum**
- **Restore natural meibum flow**
- **Better pathology management**
- **Increase glands functionality**

INTERNATIONAL GUIDE LINES ON DRY EYE (DEWS II REPORT)

International guide lines on dry eye (DEWS II Report) explains that for a good impact on the glands, heating should be done at least at 40°C. Other devices, as LipiFlow from J&J, use this concept for MGD treatment. Still the study shows how important can be meibomian gland squeezing after

heating explaining that this part is so important in order to squeeze out the melted meibum.

Activa device provides MG heating at 42°C. After few minutes it starts massaging trough pads designed to push on the lids avoiding pressing on the cornea's area and with gentle vibrations squeezing the glands out.






ACTIVA BENEFITS



PAINLESS AND RELAXING



IMMEDIATE FEELING OF RELIEF AND WELL-BEING FOR YOUR PATIENT



NO CONSUMABLES NEEDED



NO SPECIALIST REQUIRED TO OPERATE THE DEVICE



ACTIVA INCREASES THE EFFECTIVENESS OF IPL TREATMENT BY GIVING IMMEDIATE RELIEF



THE VALIDATED DISINFECTION PROCEDURES ENSURES PROTECTION FROM INFECTIONS



CUSTOMIZE NUMBER OF SESSIONS AVAILABLE FOR THE RENTAL PERIOD OF YOUR CUSTOMER



HEATING AND MASSAGE TECHNOLOGY FOR A FAST AND PRACTICAL SOLUTION

 **SBM**
Sistemi

www.sbmsistemi.com

Strada Torino, 43 - 10043 Orbassano (Torino) Italy
Tel. +39.011.19923378 - info@sbmsistemi.com