



ECO PLUS Refraction Unit

Combines modern design, comfort
and features for great value.



Control panel



Phoropter Arm



Reading lamp



Trial lens tray

- 2 Instrument sliding table top with optional electrical height adjustment.
- COMFORT Chair with arm rests and footrest, 45 degree swivel.
- Trial Lens set drawer accommodates a full trial lens set.
- Simple connection of all types of slit lamps and slit lamp illumination control.
- Convenient full featured control panel.
- Aluminium post with illumination LEDs.
- Versatile phoropter arm included.
- Clean clinical colour – off white light grey stand with mid tone grey upholstery.
- Steel base powder coated.
- Foot switch for chair height adjustment.

BOC
INSTRUMENTS

Telephone: (02) 96437888

Toll Free: 1800804331

Email: sales@bocinstruments.com.au

Website: www.bocinstruments.com.au

...Products you trust