



## **PERKINS MK3 TONOMETER**

Safety & accuracy on contact

## Perkins Mk3 The 'gold-standard'

Established for over 50 years, Perkins is synonymous with being the 'gold-standard' in hand-held applanation tonometers. The Perkins Mk3 hand-held applanation tonometer provides the same proven accuracy as other clinically-trusted Haag-Streit manufactured tonometers.

Clement Clarke originally developed the popular Perkins tonometer during the 1970s. Over time, it has been updated and improved, creating an ergonomic design, as well as giving the instrument a longer-lasting battery and enhanced internal LED illumination.

It is compatible with the Tonosafe disposable tonometer prism, which is recognised as the ideal solution for reducing the risk of cross-infection.

The Perkins Mk3 is an ideal portable applanation tonometer for use outside the traditional consulting room environment.

## Portable

The hand-held Mk3 is comfortable to hold and well-balanced. It has been specifically designed to examine patients who have often been difficult to assess by traditional applanation tonometry methods. These include; bed-ridden, wheelchair bound, paediatric and community-based patients.

## Built-in LED illumination

The innovative built-in LED illumination produces an enhanced image of the tonometer mires, permitting use regardless of ambient lighting conditions.

The LED source uses less energy and increases the operating times achieved from the integrated rechargeable battery, as well as producing lower operating temperatures, giving greater comfort to both the operator and the patient.

## Use disposable or reusable prisms

Like previous Perkins and Goldmann tonometers, the Perkins Mk3 is compatible with Tonosafe disposable prisms, so there is no need to disinfect or sterilise prisms when using outside the traditional clinic environment. In addition, the Perkins Mk3 is compatible with reusable Goldmann tonometer prisms.

The new easy-to-use 'Tonogrip' prism holder allows for straightforward placement of both sterilised reusable prisms and Tonosafe disposable prisms.

## Goldmann-compatible

The Perkins Mk3 tonometer has the unique ability to provide applanation tonometry readings that correlate strongly with the Goldmann tonometer, even with patients in a supine position.

The easy-to-use, lightweight and portable design of the Perkins Mk3 makes this tonometer an ideal choice for measurement of Intraocular Pressure (IOP).



BUILT-IN LED ILLUMINATION



EASY-TO-USE TONOGRIP PRISM HOLDER



IDEAL CHOICE FOR MEASURING IOP

## **HAAG-STREIT UK**

Edinburgh Way  
Harlow  
Essex  
CM20 2TT  
Phone (01279) 456261  
info@haag-streit-uk.com  
www.haagstreituk.com

## **Members of HAAG-STREIT Group**

HAAG-STREIT Holding AG  
www.haag-streit-holding.com

HAAG-STREIT Deutschland GmbH  
www.haag-streit.de

HAAG-STREIT AG  
www.haag-streit.com

IPro GmbH  
www.ipro.com

SPECTROS AG  
www.spectros.ch

CLEMENT CLARKE Ltd.  
www.clement-clarke.com

HAAG-STREIT Medtech AG  
www.haag-streit-medtech.com

HAAG-STREIT UK  
www.haag-streit-uk.com

HAAG-STREIT France EURL  
www.haag-streit.fr

John Weiss Ltd.  
www.johnweiss.com

HAAG-STREIT Far East  
www.haag-streit-fareast.com

HAAG-STREIT USA  
www.haag-streit-usa.com

HAAG-STREIT Surgical GmbH  
www.haag-streit-surgical.com

Reliance Medical Inc.  
www.haag-streit-usa.com

Möller-Wedel GmbH & Co KG  
www.moeller-wedel.com

Asetronics AG  
www.asetronics.ch

Möller-Wedel Optical GmbH  
www.moeller-wedel-optical.com

ComLab AG  
www.comlab.ch